

# Heart Sense

## Every Beat Counts

Hardware2 Project  
Group2 Binary



## Vision & Goal

## Features & Workflow

## Tech & Architecture

## Testing & Validation

## Conclusion & Lessons



Heart rate variability (HRV) is a **powerful indicator** of stress, recovery, and general health.



Built with Raspberry Pi Pico – **keeping cost under €10** without compromising functionality.



We believe that everyone should be able to access their physiological signals — **in real time, in a simple way.**

### Our Vision:

To make HRV tracking **accessible, affordable, and interactive** — starting with a compact device that listens to every beat and empowers every user.

## User-Friendly Design

### Single Rotary Encoder

- ✓ No extra buttons needed
- ✓ Rotate to **navigate menus**
- ✓ Press to **start/stop measuring**
- ✓ **Ideal for elderly and children**  
— intuitive, minimal, and effortless interaction




Criteria ★☆☆

## Clear Measuring & Live Wave

### Real-time interaction + data clarity

- ✓ **Live heart signal** dynamically updating
- ✓ Clear **countdown** during measurement to guide the user
- ✓ Data sent to **PC** via MQTT for deeper analysis




Criteria ★★★

## Stress Level Feedback

### Stress awareness in real time

- ✓ Uses SNS/PNS analysis **from Kubios** to estimate stress
- ✓ **Feedback** shown on OLED
- ✓ Results **stored locally** on Pico for future review



Criteria ★★★★

Vision &  
Goal

Features &  
Workflow

Tech &  
Architecture

Testing &  
Validation

Conclusion &  
Lessons

## Tired of the clunky Thonny? Magic~

✓ Local Python 3 • Thonny's Python


MicroPython (Raspberry Pi Pico) • Board in FS mode @ /dev/

MicroPython (RP2040) • Board in FS mode @ /dev/cu.usbr

Configure interpreter...

**Life saver:**

**mpremote cp \*.py :**



```
1 from piotimer import Piotimer
2 import time
3 from oled import oled
4 from basic_analysis import Basic_Analysis
5 from menu import Encoder
6 from heart_rate import CalculateHeartRate
7 from kubios import send_to_kubios, save_history
8 import ujson
9
10 rot = Encoder(10, 11, 12)
11 rot.welcome()
12 entered_first_option = False
13 hr_measurement = None
14 # for basic analysis and kubios analysis, minimum 30 seconds
15 seconds_limit = 30
16 # show countdown timer when measuring
17 timer_sec = 0
18
```

# Main & Modules

*A clean structure: main.py handles logic, other files handle the magic.*

main.py	OLED.py	Manages a single OLED instance to avoid I2C conflicts
	Menu.py	Contains all functions related to the main menu and navigation
	Heart_rate.py	Handles real-time signal acquisition and waveform plotting
	Basic_analysis.py	Performs local calculations such as Mean PPI, SDNN, etc
	MQTT.py	Manages MQTT connections and data transmission
	Kubios.py	Handles communication with Kubios via JSON upload/download

# Algorithm & Improvement (1/2)

## Smoothing the Plot Curve

- **Separate FIFO queue** stores raw data.
- Every 25 samples take average.
- Added to the main FIFO.

**Benefit:** Avoiding sampling rate issues from different sensors.

```
self.samples_unsmoothed = Fifo(25)  
self.samples = Fifo(500)
```

- **Scaled\_history list** to store the last 128 scaled samples.
- When new data arrives, the oldest data is removed.

**Benefit:** Ensuring a continuous display.

```
self.scaled_history.pop(0)  
self.scaled_history.append
```

- **Uses oled.line(x, y1, x+1, y2, 1)** to connect consecutive data points.
- Instead of oled.pixel()

**Benefit:** Eliminating the “pixelated” look.


```
for x in range(127):  
    y1 = 48 - data[x]  
    y2 = 48 - data[x+1]  
    oled.line(x, y1, x+1, y2, 1)
```

# Algorithm & Improvement (2/2)

## Improving Detection Accuracy

**Double Threshold Setting:**  
both high and low thresholds  
per second

**Benefit:** Reducing the  
occurrence of a secondary  
small peak.



```
self.threshold_high = :  
self.threshold_low = i
```

- **Data Slicing for Stability:**  
Excluding the first 4 intervals  
when send via MQTT

**Benefit:** Removing early,  
unstable data points at the  
beginning of measurement.

```
basic_analysis = Basic_Analysis()  
basic_analysis.get_result  
(hr_measurement.intervals[4:])
```

- **Heart Rate Validity Check:**  
Only calculate intervals within a  
valid BPM range (40-200 BPM)

**Benefit:** Excluding extreme values  
at the beginning to ensure the  
quality of the interval list.

```
if 40 <= bpm <= 200:  
    self.isAnalyzing = False  
    display_text = f"{bpm} BPM"  
else:  
    self.isAnalyzing = True  
    display_text = "-- analyzing"
```



Our algorithm achieved a **close match** ( $\pm 5$  BPM) with smartwatch readings under normal conditions

### 1. Comparison with Smartwatch Data

Used commercial smartwatches (e.g., Huawei, Xiaomi Watch) to compare BPM values in real time under resting conditions.

### 2. Stability Test

Tested the device across different durations (1 min, 3 min) to observe consistency in BPM readings.

### 3. Stress Testing with Noise or Motion






Simulated real-world conditions to ensure the algorithm avoids false peaks caused by small movements or signal noise.

### 4. Edge Case Detection

Tested the algorithm's response to:

- ✓ Abnormally fast/slow BPM
- ✓ Sudden signal drops
- ✓ Users with irregular waveforms (e.g., dual peaks)



Vision & Goal	Features & Workflow	Tech & Architecture	Testing & Validation	Conclusion & Lessons
 <p>Working with hardware requires patience—sensor noise and wiring issues were frequent.</p>	 <p>Importance of threshold tuning and filtering for PPG.</p>	 <p>Time synchronization between cloud and device is essential for accuracy.</p>	 <p>Visualization helps in debugging and interpretation.</p>	 <p>Task distribution and weekly planners helped keep on schedule.</p>

# Thank you!

Questions?

Hardware2 Project

Group2 Binary

Puntawat Subhamani

Olga Chitembo

Yue Zhang

